

# VERRUCA-FREEZE®

Safe, Easy & Effective Cryosurgery

www.verrucafreeze.com



Verruca-Freeze® is one of the coldest portable cryosurgical systems freezing at -70°C (-94°F). A 3-6 second spray, freezing for 20-30 seconds, makes it a quick, self-anesthetizing and easy procedure. With two forms of application: new and improved CryoCones® and CryoBuds®, Verruca-Freeze® is the clear choice for treating benign skin lesions.



## Benefits:

**Freezes at -70°C**

**1-2 minute procedure**

**Non-flammable (No Hazmat fees)**

**Self-anesthetizing**

**Minimum 5-year shelf life**

**Bloodless/Odorless**

**Portable and self-contained**

**CFC Free**

## Treatable Lesions:

warts and plantar warts	condyloma
seborrheic keratoses	granuloma annulare
skin tags (achrochordon)	prokeratosis plantaris
age spots	molluscum contagiosum
small keloids	angiomas
actinic keratoses	leukoplakia

...and more!

## Flexible Purchasing Options:

First time users should start by purchasing one of our starter kits. Replacement canisters and applicators may be purchased individually.



**CRYOSURGERY, INC.**

5829 Old Harding Pike, Nashville, TN 37205  
p. 615.354.0414 | f. 615.354.0466  
info@cryosurgeryinc.com

See back for product information and purchasing options.  
Available through medical supply distributors.

# VERRUCA-FREEZE®

Safe, Easy & Effective Cryosurgery  
www.verrucafreeze.com



## 65 Freeze Kit\*

Item #CVFK65 (162ml)

Kit includes: 162ml canister, one set of reusable CryoCones® in 5 sizes (3, 5, 7, 9 and 12mm), extender tube, cryo insulator for canister, pack of 30 disposable CryoBuds® (10 Small, 10 Medium and 10 Large), instruction manual, practice materials, training DVD & patient information sheets



## 35 Freeze Kit\*

Item #CSMK35 (95ml)

Kit includes: 95ml canister, one set of reusable CryoCones® in 5 sizes (3, 5, 7, 9 and 12mm), extender tube, cryo insulator for canister, pack of 15 disposable CryoBuds® (5 Small, 5 Medium and 5 Large), instruction manual, practice materials, training DVD & patient information sheets



## 100 Freeze Canister\*

Item #CVFL (236ml)

Replacement canister includes: 236ml canister, extender tube, instructions & patient information sheets



## 65 Freeze Canister\*

Item #CVFC65 (162ml)

Replacement canister includes: 162ml canister, extender tube, instructions & patient information sheets



## 35 Freeze Canister\*

Item #CSMC35 (95ml)

Replacement canister includes: 95ml canister, extender tube, instructions & patient information sheets



## CryoBuds®

Disposable foam tip applicators in 3 sizes (30 per pack)

Item #BUDS Small  
(for treating lesions 3mm in size)

Item #BUDM Medium  
(for treating lesions 5mm in size)

Item #BUDL Large  
(for treating lesions 9mm in size)

Item #BUDP Multi-pack  
(10 Small, 10 Medium, 10 Large)



## CryoCones®

Item #TCONES

Set of five soft-tipped, reusable cones (3, 5, 7, 9 & 12mm)



## Extender Tube

Item #TUBE

Extra tubes to insert into the end of the trigger nozzle. (pack of 5)



## Patient Information Pad

Item #PAD

Tear off information sheets that inform patients of healing process



\*Number of freezes produced by each canister will vary depending on the size of the lesion being treated and the amount of cryogen used.